


How do I order a lab on my iPad?

07/08/2024 7:33 pm EDT

In the DrChrono EHR app, you can order lab tests just as you would with the browser EHR.


To begin, navigate to a patient's chart and tap on their name. From the dropdown menu, select **Labs**.



A list of your previously ordered labs will appear. To send a new lab test, tap on the (




)button in the upper right-hand corner of the page.



A form will appear that allows you to fill out the details of your lab order. Here you will need to enter who to bill, the specimen collection date, any ICD-10 codes, notes, your office, and the lab to send to. On the requisition form,


the ordering provider will reflect the provider who is logged in. You can also select from your favorites.



You can also select the lab on the bottom half of the form. Tap on this box to search for the test you would like to conduct. You may choose multiple tests by selecting the '+' button. To save a lab as a favorite tap the star (



)icon.



After everything is filled out, select **Submit** to send your lab order. To verify that your lab was sent, you may select **Lab** from your patient chart to check that the lab you sent is now on that list.

Filters

07/27/21



Search Appointment

07/27/2021 Tue (2)



Michelle Harris
07/27 09:40AM



Jenny (Jen) Harris
07/27 12:15PM

Close Submit

Who to Bill: Send to PSC: Specimen Collection Date & Time: ★


ICD-9 + ICD-10 codes:

Note:

Ordering Provider: Office: Account:

Edit +

| | | |
|---|--|-------------|
| 1 | <input type="text" value="322000 - Comp. Metabolic Panel (14)"/> | 1 Questions |
| | <input type="text" value="Y"/> | |
| 2 | <input type="text" value="004259 - TSH"/> | 0 Questions |
| | <input type="text"/> | |
| | <input type="text"/> | |
| | <input type="text"/> | |



View All

View All

View All

Start Visit