

# OutcomeMD set up with DrChrono

Last modified on 01/13/2025 2:28 pm EST

## Setting Up DrChrono with OutcomeMD

### Download The App

Click [HERE](#) to log in via Single Sign On

- You must have Provider permission in DrChrono in order to install
- Click on **Sign in with drchrono**.

[Forgot Password?](#)



- From here, authorize the app.
- Once you've logged in and clicked on the OutcomeMD app under any patient, then you will receive a BAA to sign.
  - If you are not the individual responsible for signing BAA's for your practice, you can easily submit the proper person's email address from this screen and OutcomeMD will automatically send the BAA to this individual for signature.

---

### Login & Setup Your Center

1. After the app has been authorized, you can log into the system through DrChono.
2. To do this, log into DrChrono and get into a patient's chart by either:
  - a. **Searching for a patient** and selecting them
    - i. If on the patient list, once in, scroll to the bottom of the side left panel and look under **App Directory**.

#### App Directory


OutcomeMD
--OutcomeMD--

b. Access the patient through their appointment date within the schedule and accessing their clinical note.

ii. If within their clinical note, click on **Apps**.

3. Select **OutcomeMD**.

4. Select **Practice Overview** to be taken to OutcomeMD's portal outside of drchono.



5. Next head over to **Settings** (*Gear icon in upper right-hand corner*)

a. If you are responsible for setting up the portal, click **My Center** first.

i. Please be sure to set up the Center's:

- **Additional Settings**
  - **Enable Adverse Life Events and Patient Compliance**

▼ **Additional Settings**

- Send a SURGERY CENTER satisfaction survey with the 1st automatic postop followup
- Send a DIALYSIS CENTER satisfaction survey with the 1st automatic followup
- Send a PRIVATE PRACTICE satisfaction survey with the 1st automatic followup
- Send a question asking about ADVERSE LIFE EVENTS with each followup survey
- Send PATIENT COMPLIANCE question with each followup survey

▶ **White Labeling**

▶ **Outcome Marketing Settings**

Every patient follow up survey will include questions to assess the patients compliance with your provider's treatment plan & their adverse life events

- **White Labeling**

▼ White Labeling

Note: For White labeling purposes, please upload two logos, one 200 x 150 pixels and one 150 x 22 pixels. These will be used in patient facing areas. Accepted file formats are: PNG, JPG, JPEG, and GIF

200 x 150 pixel logo



Choose File

150 x 22 pixel logo



Choose File

▶ Outcome Marketing Settings

Cancel

Save Changes

- Add your logos (one 200 x 150 pixels and one 150 x 22 pixels).

- Outcome Marketing Settings

### PROFILE INFO

Your Specialty

My Specialty

Your Website URL

https://www.yourwebsite.com

Your Center's 'Schedule an Appointment' URL <sup>?</sup>

Enter your appointment schedule url

Email to Display (public-facing contact email shown in your marketing)

Public-facing email address

Phone Number to Display (public-facing contact phone number shown in your marketing)

Public-facing phone number

### REVIEW SITES

Patients post outcome-based reviews to your review page. Only 2 review sites may be enabled at a time.

Yelp Account Link

Choose Two

https://www.yelp.com/biz/yourlink

Google Reviews Account Link

https://business.google.com/yourlink

HealthGrades

https://healthgrades.com/yourlink

Vitals

https://vitals.com/yourlink

Enter the appropriate information.

1. Scroll to the bottom of the page and save changes

1. You can then set up the **My Profile** settings for yourself.

- a. Ask team members to set up their own **My Profile** settings.



MY PROFILE

Edit your profile

Make sure to select your **Current Center** before sending any assessments.

Current Center: Life Care

[Getting started with OutcomeMD support hub article](#)

---