

DrChrono App Unlock PIN

Last modified on 07/03/2025 6:43 pm EDT

The **DrChrono App Unlock PIN** is the four-digit code you use to access the [DrChrono EHR app](#) after a period of inactivity.



If you forget your PIN, you can log in with your username and password on the DrChrono EHR app or change your PIN from the DrChrono web application.

[View your PIN](#) | [Change your PIN](#) | [Disable the PIN on an iPad or iPhone](#) | [Enable the PIN on an iPad or iPhone](#) | [Require a PIN on an iPad or iPhone](#)

View your PIN

From the DrChrono web application


1. Select **Account** > **Account Settings** > **Profile** tab.
2. Next to **DrChrono App Unlock PIN**, select **View/Edit**.


Email Address	<input type="text" value="@drchrono.com"/>	Change Email
Password	<input type="password" value="Change Password"/>	
DrChrono App Unlock PIN	<input type="password" value="...."/>	View/Edit 4-digit numeric PIN for unlocking the mobile EHR app from inactivity.

3. Enter your DrChrono password and select **Confirm**.

From an iPad or iPhone

1. Tap **Account** > **Settings** > **DrChrono App Unlock PIN**.





2. Tap **View Your PIN**.

The PIN appears briefly.

Change your PIN

From the DrChrono web application

1. Select **Account** > **Account Settings** > **Profile** tab.
2. Next to **DrChrono App Unlock PIN**, select **View/Edit**.
3. Enter your DrChrono password and select **Confirm**.
4. Change your PIN and then select **Update Entire Profile**.

The PIN has been changed for both the DrChrono web application and the DrChrono EHR app.

From an iPad or iPhone

1. Tap **Account** > **Settings** > **DrChrono App Unlock PIN**.


2. Tap **Change Your PIN**.
3. Enter your four-digit PIN.
4. Reenter your PIN.

The PIN has been changed for both the DrChrono EHR app and the DrChrono web application.

Disable the PIN on an iPad or iPhone

1. Tap **Account > Settings > DrChrono App Unlock PIN**.
2. Tap **Disable 4-digit PIN Protection**.
3. Tap **Disable**.

The only option available under **DrChrono App Unlock PIN** is **Enable 4-digit PIN Protection**.




Enable the PIN on an iPad or iPhone

1. Tap **Account > Settings > DrChrono App Unlock PIN**.
2. Tap **Enable 4-digit PIN Protection**.
3. Enter your four-digit PIN.

4. Reenter your PIN.

The PIN has been changed for both the DrChrono EHR app and the DrChrono web application.



Require a PIN code on an iPad or iPhone

1. Tap **Account** > **Settings** > **DrChrono App Unlock PIN**.
2. Tap **Require PIN code**.
3. Select the PIN entry timeout duration before a PIN is required to log back in.

10:57 AM Thu Jul 3

82%

DashboardEHRMessagesTasksAccount

Account

< DrChrono App Unlock PIN

Settings

Practice Chat

Medical Forms Library

Help

Rate DrChrono EHR

Account Deactivation

Terms of Service

Lock

Log Out

Break the Glass

PGBlue Sky Medical

UserProvider

App Versionv11.1.3

Immediately

After 1 minute

After 5 minutes

After 15 minutes

After 1 hour

after 4 hours